


Product Data Sheet/Certificate of Analysis

HyNic-silane

Storage: Store desiccated at -20C under inert atmosphere.



C₁₈H₃₂N₄O₄Si MW 396.56

Catalog Number:	HL-1002	Unit Size:	
Lot Number:		Components:	HyNic-silane

Test	Specification	Result
Appearance (color)	Pale brownish-yellow	Passed
Appearance (form)	Solid	Passed
Purity (1H NMR)	> 90%	Passed
HPLC analysis after packaging	Conforms to structure	Passed

QC Release Date	Expiration Date:
Released By:	

Product Description

This linker should be used in conjunction with 4FB modified biomolecules for the immobilization of biomolecules to glass surfaces. The HyNic-4FB conjugation reaction is catalyzed by the addition of the TurboLink Catalyst Buffer (Cat # S-2006-105).

Comments

Store desiccated at -20C under inert atmosphere. Not for internal or external use in humans.

