



Product Data Sheet

Product: MagnaLink Amino Magnetic Beads



Storage: Store at 4°C

Catalog Number: M-1005

Diameter: 2.8 μ

Lot Number:

Expiration:

Primary amine content: XX.XX nmol/mg beads

Density: 1.9 mg/ml

Solids content: 10 mg/mL

Beads per ml: 8×10^8

Components: Uniform polymer encapsulated magnetic bead containing primary amino functional groups (-NH₂).

Size: 1.0 mL
Units: 10 mg

Product Description

A 2.8 micron uniform polymer encapsulated magnetic bead containing primary amine functional groups

Application

Amino functional groups can be activated with any water-soluble NHS-activated linker.

Comments

Not for internal or external use in humans.

