


## Product Data Sheet/Certificate of Analysis

### Sulfo-S-4FB (succinimidyl 4-formylbenzoate sodium salt)

**Storage:** Store desiccated at or below room temperature



$C_{12}H_8NO_8S.Na^+$  MW 349.25

Catalog Number:	S-1008	Unit Size:	
Lot Number:		Components:	Sulfo-S-4FB

Test	Specification	Result
Appearance (color)	White	
Appearance (form)	Solid	
Purity (1H NMR)	> 85%	
Solubility	> 50mg/mL in DMF	

QC Release Date	Expiration Date:
Released By:	

## Product Description

Sulfo-S-4FB is a water-soluble heterobifunctional linker used to convert primary amines to 4-formylbenzamide moieties.

## Comments

**Store under inert atmosphere desiccated.** Not for internal or external use in humans.

