


Product Data Sheet/Certificate of Analysis

FMOC HNA OSu (succinimidyl 6-FMOC-hydrazinonicotinate)

Storage: Store desiccated at room temperature



$C_{25}H_{20}N_4O_6$ MW 472.2

Catalog Number:	S-3001	Unit Size:	
Lot Number:		Components:	FMOC HNA OSu

Test	Specification	Result
Appearance (color)	pale yellow	
Appearance (form)	Solid	
Purity (1H NMR)	> 90%	
TLC:	Single spot	

QC Release Date	Expiration Date:
Released By:	

Product Description

6-Fmoc-HNA-Osu is a heterobifunctional linker used to convert primary amines to Fmoc-protected HyNic moieties during solid phase peptide synthesis

Comments

Not for internal or external use in humans.

