


Product Data Sheet/Certificate of Analysis

BOC HNA OSu (succinimidyl 6-BOC-hydrazino-nicotinic acid)

Storage: Store desiccated at room temperature



$C_{15}H_{18}N_4O_6$ Mol. Wt.: 350.2

Catalog Number:	S-3002	Unit Size:	
Lot Number:		Components:	BOC HNA OSu

Test	Specification	Result
Appearance (color)	pale yellow	
Appearance (form)	Solid	
Purity (1H NMR)	> 90%	
TLC:	Single spot	

QC Release Date	Expiration Date:
Released By:	

Product Description

Incorporation of BOC protected HyNic moiety on peptides or other moiety.

Comments

Not for internal or external use in humans.

