

Product Data Sheet/Certificate of Analysis

Magnetic Stand



Catalog Number:	S-6001-001		
Lot Number:	N/A	Expiration: None	
Components:	Magnetic Stand	Size: 38 mm x 24 mm x 32 mm	Units: 1

Product Description

Magnetic separation device for holding a single 1.5 mL centrifuge tube. The rack contains two embedded neodymium-iron-boron magnets used for separation of affinity magnetic beads.

Application

Used for the magnetic purification.

Comments

Warning-this product contains strong permanent magnets. People wearing a pace maker or any other magnetizable medical implant should not use this product. Implants can be affected or damaged by a strong magnetic field. Keep magnetizable tools and objects away from the working area. Keep the magnet away from credit cards and electronic recording or other electronic devices such as computers and monitors.

