



Product Data Sheet/Certificate of Analysis

HyNic-TAT-NH₂

Storage: Store desiccated at $\leq -20^{\circ}\text{C}$

HyNic-TAT.NH₂ MW 1401

Catalog Number: SP-C002

Unit Size:

Lot Number:

Components: HyNic- GRKKRRQRRRPPQ -NH₂

Test	Specification	Result
Purity (LC)	>90%	
Uv/Vis	214nm	
Mass (positive scan)	1852 (peptide), 1892 (peptide+acetone)	
Appearance- Color	White	
Appearance- Form	Solid	

QC Release Date

Expiration

Date:

Released By:

Product Description

Pre-activated Peplink cell penetrating peptide used for intracellular delivery

Comments

Store desiccated at $\leq -20^{\circ}\text{C}$. Not for internal or external use in humans.

The products offered here are for research use only. Any commercial application will require a license from Solulink. The Solulink Conjugation System is patented and has additional patents pending. Please contact Solulink for information regarding licensing information. No license is granted or implied to any patents to technologies for which the end user applies our products.

