

Product Data Sheet

HyNic/PEG2/Hexa-His

Compound: SP-E001 (HyNic/PEG2/Hexa-His)

Storage: Store desiccated at < -20° C

HyNic/PEG2/HHHHHH-NH2 MW 1133.47– white solid

Catalog Number:	SP-E001	Lot Number:	
Components:	HyNic/PEG2/Hexa-His	Unit Size:	

Test	Specification	Result
Purity (LC)	>90%	
Uv/Vis	214nm	
Mass (positive scan)	1132.84 (peptide), 1174.84 (peptide+acetone)	
Appearance- Color	Pale yellow	
Appearance- Form	Solid	

QC Release Date		Expiration Date:
Released By:		

Application

The linkable HyNic -(His)₆ epitope tag provides a non-biotin based, affinity tagging method for the soft-release or detection of tagged non-recombinant proteins. SoluLink is the first to provide pre-activated, linkable peptides for tagging native, recombinant proteins.

Comments

Store desiccated < -20°C. Not for internal or external use in humans.

