



Product Data Sheet

S-Tag-HyNic

Compound: SP-E002 (HyNic/PEG2/S-Tag)

Storage: Store desiccated at $< -20^{\circ}\text{C}$

(HyNic) (PEG2) SDMHQREFKAAAYGK-NH2

MW1133.47- white solid

Catalog Number:	SP-E002-010	Purity:	$> 90\%$ by ^1H NMR
Lot Number:		Expiration:	
Components:	HyNic/PEG2/S-Tag	Unit Size:	1.0 mg

Application

The linkable HyNic -S- epitope tag provides a non-biotin based, affinity tagging method for the soft-release or detection of tagged non-recombinant proteins. Solulink is the first to provide pre-activated, linkable peptides for tagging native, recombinant proteins, oligonucleotides, or other biological molecules.

Comments

Store under inert atmosphere desiccated $< -20^{\circ}\text{C}$. Not for internal or external use in humans.

